

Diversey



MoonBeam™3 Disinfection Technology

On-Demand Targeted Disinfection System

MoonBeam™3 delivers powerful UV-C light to provide fast, on-demand broad area disinfection of high-touch surfaces in as little as 3 minutes.

Features & Benefits

- Three individually adjustable articulating heads, each generating a broad umbrella of intense UV-C light. Cost-effective operation – low operating cost and minimal maintenance
- Base unit - locking casters allow easy movement and positioning anywhere. UV-C heads unfold for positioning in seconds
- Remote cover - controls operation and acts as a sensor during disinfection to ensure safety. Cover also protects light heads and enables convenient transport from room to room. Optional multi-access point motion sensor also available
- Added Assurance - provides additional disinfection for surfaces that may have been missed in cleaning
- Ideal solution for frequently-touched, hard-to-clean surfaces and intricate, hard-to-wipe equipment. Fume-free, chemical-free disinfection

Applications

- Suited for whole-room or targeted disinfection of high-touch surfaces and equipment including: patient bed area, bathrooms, fixtures, work stations on wheels (WOWs), etc.





MoonBeam™3 Disinfection Technology

On-Demand Targeted Disinfection System




Product Specifications

Description	MoonBeam™3
UV-C vertical dosing height range	6" (15 cm) to 84" (210 cm) high
Cycle time in seconds	90, 180, 300, 600
Device weight	39 pounds (17 kg) complete 27 pounds (13 kg) Base Unit
Environmental	No ozone generation, high output UV-C lamps
Device Size	15" x 44" (38 cm x 112 cm)
Acoustic Noise	Negligible
Power	100-120V~, 3.2A, 50/60Hz 200-240V~, 1.6A, 50/60Hz

Replacement Parts (Available from Plan B, Inc.)

Product Code	Description
SP-00695	MoonBeam™3 Remote Cover
SP-00696	MoonBeam™3 UV-C Head Reflector
SP-00697	MoonBeam™3 UV-C Head Endcap
SP-00700	MoonBeam™3 Fuse
SP-00720	MoonBeam™3 Cable, Line Cord, North America
SP-00701	MoonBeam™3 Packaging

MoonBeam™3 Device Components

Product Code	Description	Details
MOONBEAM3	MoonBeam™3 UV-C Disinfection Technology Device	MoonBeam™3 is designed for fast, on-demand disinfection of high-touch surfaces 
MOON3-SENTRY	MoonBeam™3 Sentry Assist Technology	Additional 360-degree motion sensor remotely terminates MoonBeam™3 cycle if human movement is detected in disinfection zone 
MOON3-UVCBULBS	MoonBeam™3 Bulb Set (3 bulbs)	3-UV-C Replacement Bulbs for MoonBeam™3 

Safe handling

Please make sure your employees read and understand the product label and Safety Data Sheet before using this product. The label contains directions for use; and both the label and SDS contain hazard warnings, precautionary statements and first aid procedures. SDS are available online at www.diversey.com or by calling 888.352.2249. Improper use or dilution may result in damage to surfaces and may result in health and physical hazards that match those of the concentrate.

Ask about our Return to Factory Warranties available for up to 5 years.
Contact your Diversey Sales Representative.

www.diversey.com

© 2018 Diversey, Inc. All Rights Reserved. PSS90612 en (19/101)