

# 3Gen COMPARISON CHART

	<a href="#">request a visit</a>	<a href="#">click to order</a>	<a href="#">click to order</a>	<a href="#">click to order</a>	<a href="#">click to order</a>	<a href="#">click to order</a>	<a href="#">click to order</a>	<a href="#">click to order</a>	<a href="#">click to order</a>	<a href="#">click to order</a>	<a href="#">click to order</a>	<a href="#">click to order</a>
<b>Click to order&gt;&gt;</b>												
<b>Ask about volume discounts!</b>	<b>Lumio 108-0398</b>	<b>Lumio UV 218-0117</b>	<b>Lumio S 120-3518</b>	<b>DL3 711-0040</b>	<b>DL3N 121-2527</b>	<b>DL1 117-4300</b>	<b>DL100 711-0188</b>	<b>Carbon 110-6516</b>	<b>DL II Pro 711-4811</b>	<b>DL II Pro HR 711-4286</b>	<b>2 Pro Hybrid 108-7376</b>	<b>DL4 711-0117</b>
<b>LEDs</b>	40 white	40 UV LEDs	24P + 8NP white	21P + 7NP white	12P + 6NP, 10 PigmentBoost™	4 white	8 white	8P + 8NP white	16P + 16 NP white	32 white	24P + 8NP white	12 P + 6 NP, 6 Pigment Boost Plus
<b>Polarization System</b>	Permanent Cross-Polarization	None	Cross- and No Polarization	Cross- and No Polarization	Cross- and No Polarization	Cross- and No Polarization	Permanent Cross-Polarization	Cross- and Linear Polarization	Cross- and Linear Polarization	Permanent Cross-Polarization	Cross- and No Polarization	Cross- and No Polarization
<b>Lens   Optics</b>	75 mm, 2x	75 mm, 2x	40mm, 4x	Four-element, 25mm, 10x	Four-element, 25mm, 10x	15mm, 10x	15mm, 10x	15mm, 10x	Three-element, 25mm, 10x	Three-element, 25mm, 10x	Three-element, 25mm, 10x	Four-element, 30mm, 10x
<b>Power Source</b>	4 Lithium AA	4 Lithium AA	Rechargeable USB Cable	Rechargeable Lithium	Rechargeable Lithium	Rechargeable Lithium	2CR5 Lithium Battery	2CR5 Lithium Battery	Rechargeable Lithium	Rechargeable Lithium	Rechargeable Lithium	Rechargeable Lithium
<b>Camera Adaptable</b>	No	No	No	Yes  Nikon1  iPhone 5S/5C   Galaxy S3/S4   iPads	Yes  Nikon1  iPhone 5S/5C   Galaxy S3/S4   iPads	Yes   all iPhone   iPads   Galaxy S3/S4	No	No	Yes   all iPhone   iPads   Galaxy S3/S4	Yes   all iPhone   iPads   Galaxy S3/S4	Yes   all iPhone   iPads   Galaxy S3/S4	Yes  Nikon1  all iPhone   Galaxy S3-S6   iPads
<b>No-Oil Lesion Imaging</b>	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
<b>Immersion Fluid Imaging</b>	No	No	No	Yes	Yes	Yes	No	No	No	No	Yes	No
<b>Superficial Skin Imaging</b>	No	No	No	No	No	No	No	Yes	Yes	No	No	No
<b>Color-Based Skin Imaging</b>	No	No	No	No	No	No	No	No	No	No	No	Yes
<b>Measuring Scale (mm)</b>	No	No	No	10 mm	10 mm	10 mm	No	No	10 mm	10 mm	10 mm	10 mm
<b>Body Style</b>	Polycarbonate white	Polycarbonate white	Gorilla Glass	Aluminum, black & silver	Aluminum, black & silver	Aluminum, black	Polycarbonate, white	Polycarbonate, black	Polycarbonate, black	Polycarbonate, black	Polycarbonate, navy	Aluminum, black & silver
<b>Carrying Feature</b>	Neoprene Pouch	Neoprene Pouch		Black Leather Case   Lanyard	Black Leather Case   Lanyard	Pocket Clip	<i>Optional</i>	Silicone Sleeve   Lanyard	Silicone Sleeve   Lanyard	Silicone Sleeve   Lanyard	Silicone Sleeve   Lanyard	Black Leather Case   Silicone Sleeve   Lanyard