

seca®

HENRY SCHEIN®
MEDICAL

LIMITED TIME!
**MANUFACTURER
REBATE**

See details inside...

SCALING UP

**WEIGHING &
MEASURING REDEFINED**

seca SCALE-UP LINE



Play video to see the scale in action!

Built for Durability Designed for Safety Ready for High Traffic

MOBILITY MEETS DURABILITY

The optimized mobile and easy-to-move model of the seca Scale-up Line incorporates a digital height rod and optional rechargeable battery pack. Digital height measuring results appear in the display automatically with millimeter precision.

The seca Scale-up Line is purpose-built to endure the rigors of daily medical work with resilience. Its heavy-duty, robust design ensures consistent accuracy even after years of relentless service, providing healthcare teams with trusted reliability.

Large roller castors ensure easy mobility for nursing professionals on the go.



The optional rechargeable battery pack provides more than 8 hours of continuous weighing and measuring.



Obtain weight and digital height conveniently at each patient room, bedside, or wherever needed.



The rock-solid steel handrail enhances stability to reduce fall risk.




DESIGNED TO THRIVE

in Demanding Healthcare Environments

- ① Ultrasound height measuring
- ② Sturdy aluminum column
- ③ Customizable touch-display
- ③ Voice-guided automation
- ④ Robust steel handrail
- ⑤ Patient and user identification*
- ⑥ Optional Rechargabel battery pack
- ⑦ Heavy-duty wheels
- ⑧ Break-resistant Bearclaw glass
- ⑧ Easy to clean surfaces

- + Capacity of 800lb / 360kg
- + Fine 0.1oz / 50 g graduation

* Optional barcode scanner mount seca 463, and barcode scanner from a third-party vendor.




The seca Scale-up Line redefines medical scales' standards, delivering precision and versatility that will reshape your expectations. Its innovative features and unmatched durability set a new benchmark in the industry, making it the go-to choice for healthcare professionals.


One Scale to Master all Weighing Challenges

The seca Scale-up Line is your versatile solution, perfect for a wide range of patients and healthcare scenarios. Whether you're working with pediatric or geriatric patients, need mobile functionality, or require seamless EMR integration, this scale was designed to excel in every environment.


Quick digital height and weight measurement during triage.



The optional panda motif makes taking children's measurements lots of fun.



The ultrasonic measuring station model can be tailored for young patients with a fun animal sticker accessory.



Bedside support for convenient weighing.

Seamless Direct EMR-validated Integration



Enhance efficiency, eliminate transcription errors, and foster a zero-incident environment with the seca Scale-up Line. Seamlessly integrate patient measurements into your EMR* System. Imagine the relief for staff as patient weight and height data is automatically captured. Contact us for more information on our integration solutions.

→ seca.com/integration



* EMR integration features are functional exclusively when the device is connected to the AC Adapter and are not applicable when powered by a mobile battery pack.



The seca ID-Display shows patient and user IDs and transfers measurement data with the push of a button.



Optional barcode scanner mount seca 463, and barcode scanner from a third-party vendor.

TAILOR YOUR SOLUTION

Find the perfect Scale-up Line model for your needs



For more information, please contact your Henry Schein Representative or call Henry Schein at 1-800-P-SCHEIN (1-800-772-4346).

seca north america | west
13601 Benson Avenue Chino,
CA 91710
USA

T +1800-542-7322
F +1888-705-7397
E info.us@seca.com

seca.com/scaleup

LIMITED TIME!
\$200 Manufacturer
Rebate!*



Features	Handrail scale with ID-Display™	Mobile handrail scale with digital measuring rod and ID-Display™	Ultrasonic measuring station with ID-Display™	Ultrasonic measuring station with ID-Display™ and handrail
+ Flat scale seca 654 Capacity: 800 lbs / 360 kg or 360 k Graduation: 0.1 lb / 50 g or 50 g	●	●	●	●
+ Handrail seca 455	●	●		●
+ Digital height rod seca 254: 4 – 90 inch / 10 – 228 cm (digital measuring: 48 – 90 inch / 122 – 228 cm)		●		
+ Ultrasonic height measurement seca 257: without handrail: 24 – 87 inch / 60 – 220 cm with handrail: 40 – 87 inch / 100 – 220 cm			●	●
OPTIONAL ACCESSORIES + Rechargeable battery pack seca 453 Can be retrofitted to Scale-up Line devices with serial number 10000001350445 or higher	●	●		
+ Animal sticker seca 487			●	●
+ Panda seca 459 Not in combination with battery pack	●	●		
+ EMR integration The features are functional exclusively when the device is connected to the AC Adapter and are not applicable when powered by a mobile battery pack	●	●	●	●
Order code	HS Item # 138-3144	HS Item # 145-9845	ON SMM IUT NN	ON SHM IUT NN

*Manufacturer rebate valid for mobile handrail scale with digital measuring rod (HS Item # 145-9845) sold within the specified promotional period (July 1, 2024 - September 30, 2024) by Henry Schein Medical in the US. This offer may not be combined with any other offers or promotions. To redeem your manufacturer rebate, contact seca at info.us@seca.com

